Diagnostic yield, complications, and downstream lung cancer treatment after shape-sensing robotic-assisted bronchoscopy (ssRAB) compared to transthoracic needle aspiration biopsy

INTUITIVE

Sam E. Wing¹, Adam Gaertner¹, I-Fan Shih¹, Yuki Liu¹ ¹Intuitive Surgical, Sunnyvale, CA USA



BACKGROUND

For patients with suspected lung cancer nodules, shapesensing robotic-assisted bronchoscopy (ssRAB) may offer improved diagnostic yields with a strong safety profile.^{1,2} This could improve access to curative treatments for patients despite the higher upfront cost of a navigational bronchoscopic biopsy compared to percutaneous biopsy methods. Faster access to treatment can also prevent adverse outcomes like upstaging and cancer progression, which can negatively impact patient survival. 3

Using a large, US-based RWE database, we analyze the biopsy results and downstream lung cancer treatments of patients undergoing either transthoracic needle aspiration (TTNA) or ssRAB lung nodule biopsies.



Data Source:

Premier Healthcare Database

Patient-level data from hospital discharges in the United States

Patient Identification:

TTNA Patients: CPT & ICD 10-PCS codes* ssRAB: Billing text strings (see case identification) Solitary pulmonary nodules: ICD 10-CM codes*

Follow Up:

Patients were followed for 1 year after their index biopsy encounter

Statistical Analysis:

Propensity score matching was performed to compare complication rates. The propensity score was built using:

- Admission type
- Age
- Gender Race/ethnicity
- Charlson Comorbidity Index
- Obesity
- Physician volumes

STRENGTHS + LIMITATIONS

Strengths:


 Identification of ssRAB patients in a large, RWD source is a novel approach & helps generalize findings from smaller clinical studies


Limitations:

- Claims database has limited clinical information, including nodule size and tumor pathology
- Lost-to-follow-up occurred if patients received treatment at an outside facility

ssRAB CASE IDENTIFICATION

ssRAB cases were identified using the following text strings found in bronchoscopy billing files:





RESULTS

Patient Population:

Among 159,961 patients receiving a biopsy, **2.8% (n=4,337) underwent an ssRAB** and 97.2% (n=152,594) underwent TTNA.

Diagnostic Results:

62% of ssRAB patients vs 52% of TTNA patients had a diagnostic result with a malignant pathology. 38% of ssRAB malignancies and 27% of TTNA malignancies received cancer treatment in the database.

Downstream Cancer Treatments:

Surgery (27.9%) or radiotherapy (27.4%) were the most common treatment approaches, with surgery and radiotherapy occurring 2.3x and 1.4x more often among ssRAB patients vs. TTNA patients.

Complications:

The unadjusted odds ratio of a pneumothorax was 9.13 (p < 0.001) comparing TTNA to ssRAB outpatient cases. The unadjusted odds ratio of a pneumothorax with a chest tube intervention was 13.07 (p < 0.001) comparing TTNA to ssRAB outpatient cases. This relationship was materially unchanged after propensity score matching.

Pneumothorax Diagnosis (2020 - 2023)		
	lon	TTNA
Total Procedures w/ SPN Diagnosis (n)	3,823	61,123
Inpatient (n, %)	19 (0.5%)	14,886 (24.4%)
Outpatient (n, %)	3,804 (99.5%)	46,237 (75.6%)
Pneumothorax diagnoses, 0-3 days post bx (n, %)	37 (1.0%)	4,112 (8.9%)
7-day reencounter rate (n, %)	2 (5.4%)	1,439 (35.0%)
Reencounter LOS (Mean, SD)	0 (0)	0.15 (0.75)
Pneumo. dx w/ chest tube, 0-3 days post bx (n, %)	10 (0.3%)	1,589 (3.4%)
7-day reencounter rate (n, %)	2 (20.0%)	576 (36.2%)
Reencounter LOS (Mean, SD)	0 (0)	0.13 (0.63)

CONCLUSION

- Patients undergoing an ssRAB were more likely to receive a malignant diagnosis, possibly improving their ability to expedite entry into the treatment path and avoid the negative consequences of delaying cancer care
- ssRAB was also linked to a higher continuity of care, where patients were more likely to receive treatment within the same hospital
- TTNA encounters were significantly more likely to result in a complication, demonstrating ssRAB's strong safety profile
- Future studies should evaluate the impact of ssRAB on other adverse biopsy outcomes (e.g hemorrhage and repeat biopsies), especially among patients with small lung nodules



0BD34ZX, 0BD44ZX, 0BD54ZX, 0BD64ZX, 0BD74ZX, 0BD84ZX, 0BD94ZX, 0BDB4ZX, 0BDC4ZX, 0BDD4ZX, 0BDF4ZX, 0BDH4ZX, 0BDJ4ZX, 0BDL4ZX, 0BDL4ZX, 0BJ04ZZ, 0BJK4ZZ, 0BJK4ZZ, 0BJL4ZZ, 0W993ZX, 0W994ZX, 0W9B3ZX, 0W9B4ZX, 0W9BAXX, 0W9 Pneumothorax Diagnosis Codes: J93, J93.0, J93.1, J93.11, J93.12, J93.8, J93.82, J93.83, J93.9, J95.811, J95.812,

Yarmus L, Akulian J, Wahidi M, et al. A Prospective Randomized Comparative Study of Three Guided Bronchoscopic Approaches for Investigating Pulmonary Nodules: The PRECISION-1 Study. Chest. 2020;157(3):694-701. Simoff MJ, Pritchett MA, Reisenauer JS, et al. Shape-sensing robotic-assisted bronchoscopy for pulmonary nodules: initial multicenter experience using the Ion™ Endoluminal System. BMC Pulmonary Medicine. 2021;21(1):322. Mohammed N, Kestin LL, Grills IS, et al. Rapid Disease Progression With Delay in Treatment of Non-Small-Cell Lung Cancer. International Journal of Radiation Oncology*Biology*Physics. 2011;79(2):466-472.